Case: Normal Operation (Beginning a Session) MainWindow sessionTimer breathTimer This diagram shows the normal operation of the HeartWave device. The User is in the main menu where the first option is to begin a session. The session ends when the select button is pressed. After, the log of the session is displayed on the screen. press() handleSelectButtonPress() startSession() setBreathInterval() setChallengeLevel() startTimer(5000) getCoherenceScore() ordingLengthofSession() getFiveLatestPlotPoints() press() Log() creates Log object getLog(logNumber) returns Log object getPlotPoints() Log of session is displayed on screen