## Case: Change Breath Interval and Challenge Level Settings downButton selectButton MainWindow Device User This diagram shows how a user can use the UI buttons to change the breath interval and challenge level settings. The user is on the main menu, will go into the settings menu, and then select the option they want to change. press() handleSelectButtonPress() updateMenuList(settings) updates menu UI to show settings [case 1: Challenge Level Change] the challenge level option is the first on the settings screen press() handleSelectButtonPress() updates menu UI to show challenge level options updateMenuList(challengelevel) loop [until user gets to desired challenge level option] press() handleDownButtonPress() press() handleSelectButtonPress() setChallengeLevel() [case 2: Breath Interval Change] the breath pacer option is the second on the settings screen press() press() handleSelectButtonPress() updateMenuList(breathpacer) updates menu UI to show breath pacer options [until user gets to desired breath interval option] loop press() handleDownButtonPress() press() handleSelectButtonPress() setBreathPace()